

© +31(0)297 286848 +31(0)297 287570

info@biotrading.com www.biotrading.com

K.v.K. nr. 30059547

Medium	DG-18 AGAR
IVICUIUIII	DG-10 AGAR

Recipe number 103

Dichloran 18% glycerol agar Also known as

Enumeration of yeasts and moulds in food and animal feeding stuff with a water activity ≤ Intended use

Typical formula (g/l)* Peptone 5.0; Glucose 10.0; Potassium dihydrogen phosphate 1.0; Magnesium sulfate

0.5; Dichloran 2 mg; Agar 15.0; Chloramphenicol 0.1; Glycerol 220.0

ISO 21527-2 Formula conform

Interpretation Count the number of yeasts and/or moulds present and consider the dilution factor in

determining the yeast and/or mould counts per gram or millilitre of sample being tested.

Quality Control

Sterility test 7 days at 25°C

Final pH at 20-25°C 5.6 ± 0.2

Appearance Beige, trace to slightly hazy gel

Growth performance

Control strains	Test method	<u>Conform</u>	<u>Incubation</u>	<u>Criteria</u>
S. cerevisiae ATCC 9763	80-120 cfu	ISO 11133	5d 25°C	PR ≥ 0.5 (ref. SDA), characteristic colonies
A. caesiellus ATCC 42693	80-120 cfu	ISO 11133	5d 25°C	PR ≥ 0.5 (ref. SDA), characteristic colonies
W. mellicola ATCC 42694	80-120 cfu	ISO 11133	5d 25°C	PR ≥ 0.5 (ref. SDA), characteristic colonies
E. coli ATCC 25922	≥10⁴ cfu	ISO 11133	5d 25°C	no growth
A. brasiliensis ATCC 16404	streak		2d 37°C	growth
S. cerevisiae ATCC 4098	streak		2d 37°C	growth
E. coli ATCC 25922	streak		1d 37°C	total inhibition

Remarks

n.a.

Storage instructions : petridishes 2-15°C; petridishes with extended shelf life (-MF, -TF, -QC, -TR) 12-25°C; tubes/bottles/jars/bags 2-25°C;

chromogenic products P090KP/R055PV: 2-8°C; minimize exposure to light.

Version 13, 17-07-2020



^{*} Adjusted and/or supplemented as required to meet performance criteria